


Anti-Von Willebrand Factor (Vwf) Rabbit Polyclonal Antibody For Wb/Ihc-P/Elisa - P04275

Item Number: CM0013555



Introduction

Rabbit polyclonal antibody targeting human von Willebrand factor (VWF), immunogen amino acids 2645-2813. Validated for WB, IHC-P, and ELISA. Suitable for hemostasis research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	von Willebrand factor (VWF) Rabbit pAb
Remarks/Alias	VWD; F8VWF; von Willebrand factor (VWF)
Species	Human
GeneID (Human)	7450
GeneID	7450
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2645-2813 of human von Willebrand factor (VWF) (NP_000543.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P04275
Protein Weight	309kDa
Shipping	Ice bag