

Anti-Vprbp Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y4B6

Item Number: CM0020493



Introduction

Rabbit polyclonal antibody against human VPRBP (DCAF1), a substrate receptor of CUL4A-RBX1-DDB1 E3 ligase. Validated for WB and ELISA, with mouse cross-reactivity. Suitable for cell cycle, epigenetics, and host-virus interaction studies.

[Learn More](#)

Parameter	Value
Product Name	VPRBP Rabbit pAb
Remarks/Alias	RIP; VPRBP
Species	Human
Gene ID	9730
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1058-1396 of human VPRBP (NP_055518.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9Y4B6
Protein Weight	118kDa/168kDa/169kDa
Shipping	Ice bag