

Anti-Vps11 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9H270

Item Number: CM0003904

Introduction

VPS11 Rabbit monoclonal antibody, validated for Western blot, IHC-P, and ELISA applications. Detects human and mouse VPS11, a core component of HOPS/CORVET complexes in endosomal trafficking and autophagy. Ideal for IVD research and assay development.

[Learn More](#)

Parameter	Value
Product Name	VPS11 Rabbit mAb
Remarks/Alias	END1, PEP5, HLD12, RNF108, hVPS11, DYT32, VPS11
Species	Human
Gene ID (Human)	55823
Gene ID	55823
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human VPS11 (Q9H270).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H270
Protein Weight	108kDa
Shipping	Ice bag