

Anti-Vps13A Polyclonal Antibody For Wb And Elisa - Q96RI7

Item Number: CM0019986

Introduction

Rabbit polyclonal antibody targeting human VPS13A (chorein), a lipid transfer protein at organelle contact sites. Validated for WB and ELISA in human samples, with a molecular weight of 360 kDa.

[Learn More](#)

Parameter	Value
Product Name	VPS13A Rabbit pAb
Remarks/Alias	CHAC; BLTP5A; CHOREIN; VPS13A
Species	Human
GeneID	23230 (Human)
Immunogen	Recombinant protein corresponding to amino acids 350-600 of human VPS13A (NP_150648.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96RL7
Protein Weight	360kDa
Shipping	Ice bag