

Anti-Vps16 Rabbit Polyclonal Antibody For Wb And Elisa - Q9H269

Item Number: CM0020319



Introduction

Rabbit polyclonal antibody against VPS16, validated for WB and ELISA. Cross-reacts with human, mouse, rat. Ideal for studying endosomal trafficking, autophagy, and lysosomal disorders.

[Learn More](#)

Parameter	Value
Product Name	VPS16 Rabbit pAb
Remarks/Alias	DYT30; hVPS16; VPS16
Species	Human
GeneID (Human)	64601
GeneID	64601
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 540-839 of human VPS16 (NP_072097.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H269
Protein Weight	95kDa
Shipping	Ice bag