

Anti-Vps18 Rabbit Monoclonal Antibody For Wb And Elisa - Q9P253

Item Number: CM0003210

Introduction

High-quality rabbit monoclonal antibody targeting human VPS18, validated for Western blot and ELISA. Suitable for research on endosomal trafficking, autophagy, and lysosomal function. Species reactivity: human.

[Learn More](#)

Parameter	Value
Product Name	VPS18 Rabbit mAb
Remarks/Alias	PEP3; VPS18
Species	Human
GeneID (Human)	57617
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human VPS18 (Q9P253)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9P253
Protein Weight	110 kDa
Shipping	Ice bag