

Anti-Vps18 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9P253

Item Number: CM0023508

Introduction

VPS18 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets vacuolar protein sorting-associated protein 18, a HOPS/CORVET component in lysosomal trafficking. Ideal for autophagy and endocytosis research.

[Learn More](#)

Parameter	Value
Product Name	VPS18 Rabbit pAb
Remarks/Alias	PEP3; VPS18
Species	Human
Gene ID (Human)	57617
Gene ID	57617
Immunogen	Synthetic Peptide A synthetic peptide of human VPS18
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P253
Protein Weight	110kDa
Shipping	Ice bag