


Anti-Vps24 Monoclonal Antibody For Wb - Q9Y3E7

Item Number: CM0002997



Introduction

Rabbit monoclonal antibody against human VPS24 (CHMP3), validated for WB and ELISA. Detects human, mouse, and rat VPS24 with a predicted molecular weight of 25kDa.

[Learn More](#)

Parameter	Value
Product Name	VPS24 Rabbit mAb
Remarks/Alias	NEDF; VPS24; CGI-149
Species	Human
GeneID (Human)	51652
GeneID	51652
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 65-220 of human VPS24 (Q9Y3E7).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y3E7
Protein Weight	25kDa
Shipping	Ice bag