


Anti-Vps25 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Brg1

Item Number: CM0018127



Introduction

Rabbit polyclonal antibody against human VPS25 (ESCRT-II subunit), validated in WB, IF/ICC, ELISA, with cross-reactivity to mouse and rat. Ideal for studying endosomal sorting and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	VPS25 Rabbit pAb
Remarks/Alias	DERP9; EAP20; FAP20; VPS25
Species	Human
GeneID (Human)	84313
GeneID	84313
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-176 of human VPS25 (NP_115729.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BRG1
Protein Weight	21kDa
Shipping	Ice bag