


Anti-Vps26A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - 075436

Item Number: CM0015666



Introduction

VPS26A rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Recognizes human, mouse, and rat VPS26A, a key retromer component for endosomal cargo sorting and recycling.

[Learn More](#)

Parameter	Value
Product Name	VPS26A Rabbit pAb
Remarks/Alias	HB58; PEP8A; VPS26; Hbeta58; VPS26A
Species	Human
Gene ID (Human)	9559
Gene ID	9559
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 256-327 of human VPS26A (NP_004887.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075436
Protein Weight	38kDa
Shipping	Ice bag