


Anti-Vps28 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Uk41

Item Number: CM0027422



Introduction

Rabbit polyclonal antibody targeting human VPS28, validated for Western blot, immunofluorescence/immunocytochemistry, and ELISA. Cross-reacts with mouse and rat. Ideal for ESCRT-I complex and vesicular trafficking research.

[Learn More](#)

Parameter	Details
Product Name	VPS28 Rabbit pAb
Remarks/Alias	VPS28
Species	Human
GeneID (Human)	51160
GeneID	51160
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-221 of human VPS28 (NP_057292.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UK41
Protein Weight	25kDa
Shipping	Ice bag