


Anti-Vps28 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Uk41

Item Number: CM0003742



Introduction

Rabbit monoclonal anti-VPS28 antibody validated for WB, IHC-P, ELISA. Recognizes human, mouse, rat VPS28. Ideal for studying ESCRT-I complex in endosomal trafficking.

[Learn More](#)

Parameter	Value
Product Name	VPS28 Rabbit mAb
Remarks/Alias	VPS28; ESCRT-I subunit
Species	Human
GeneID (human)	51160
GeneID	51160
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 122-221 of human VPS28 (Q9UK41).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UK41
Protein Weight	25kDa
Shipping	Ice bag