


Anti-Vps29 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Ubq0

Item Number: CM0015670



Introduction

Rabbit polyclonal antibody targeting human VPS29 (SWISS: Q9UBQ0) for WB, IF/ICC, ELISA applications. Validated in human, mouse, rat. Involved in endosomal recycling and retromer complex. Ideal for studying retromer and retriever complexes, cargo sorting, and transcytosis.

[Learn More](#)

Parameter	Value
Product Name	VPS29 Rabbit pAb
Remarks/Alias	DC7; DC15; PEP11; VPS29
Species	Human
GeneID (Human)	51699
GeneID	51699
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-186 of human VPS29 (NP_476528.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBQ0
Protein Weight	21kDa
Shipping	Ice bag