


Anti-Vps35 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Qk1

Item Number: CM0027992



Introduction

Anti-VPS35 Rabbit Polyclonal Antibody validated for WB and ELISA. Targets human and mouse VPS35 (92 kDa, UniProt Q96QK1). A core retromer complex subunit essential for endosomal cargo sorting and linked to Parkinson disease. Suitable for intracellular trafficking research.

[Learn More](#)

Parameter	Value
Product Name	VPS35 Rabbit pAb
Remarks/Alias	MEM3, PARK17, VPS35
Species	Human
Gene ID	55737
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human VPS35 (NP_060676.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96QK1
Protein Weight	92kDa
Shipping	Ice bag