


Anti-Vps37B Polyclonal Antibody For Wb, Elisa - Q9H9H4

Item Number: CM0018888



Introduction

High-quality rabbit polyclonal antibody targeting human VPS37B, validated for WB and ELISA. Reacts with human, mouse, and rat samples. Ideal for ESCRT-I complex and endosomal sorting research.

[Learn More](#)

Parameter	Specification
Product Name	VPS37B Rabbit pAb
Remarks/Alias	VPS37B
Species	Human
GeneID (human)	79720
GeneID	79720
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-170 of human VPS37B (NP_078943.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H9H4
Protein Weight	31kDa
Shipping	Ice bag