


Anti-Vps39 Polyclonal Antibody For Wb, If/Icc, Elisa - Q96Jc1

Item Number: CM0018356



Introduction

High-quality rabbit polyclonal antibody targeting human VPS39. Validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for studies on vesicle trafficking, autophagy, and TGF-beta signaling.

[Learn More](#)

Parameter	Specification
Product Name	VPS39 Rabbit pAb
Remarks/Alias	TLP; VAM6; hVam6p; VPS39
Species	Human
GeneID (Human)	23339
GeneID	23339
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 606-875 of human VPS39 (NP_056104.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96JC1
Protein Weight	102kDa
Shipping	Ice bag