

Anti-Vps4A Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Un37

Item Number: CM0025188



Introduction

VPS4A Rabbit Polyclonal Antibody (pAb) validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat VPS4A. Ideal for endosomal trafficking, cytokinesis, and viral budding research.

[Learn More](#)

Parameter	Value
Product Name	VPS4A Rabbit pAb
Remarks/Alias	SKD1, SKD2, VPS4, SKD1A, CIMDAG, VPS4-1, VPS4A
Species	Human
GeneID (human)	27183
GeneID	27183
Immunogen	Synthetic peptide / Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human VPS4A (NP_037377.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UN37
Protein Weight	49kDa
Shipping	Ice bag