

Anti-Vps54 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9P1Q0

Item Number: CM0016337

Introduction

Rabbit anti-VPS54 polyclonal antibody validated for WB, IF/ICC, and ELISA. Detects human and mouse VPS54. Suitable for studying retrograde transport, GARP complex function, and tumor antigen research.

[Learn More](#)

Parameter	Value
Product Name	VPS54 Rabbit pAb
Remarks/Alias	WR; HCC8; SLP-8p; VPS54L; hVps54L; PPP1R164; VPS54
Species	Human
GeneID (Human)	51542
GeneID	51542
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 420-720 of human VPS54 (NP_001005739.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9P1Q0
Protein Weight	111kDa
Shipping	Ice bag