

Anti-Vsig1 Rabbit Polyclonal Antibody For Wb, Elisa - Q86Xk7

Item Number: CM0014210

Introduction

Rabbit polyclonal antibody targeting human VSIG1/GPA34, validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for cancer and gastrointestinal tissue studies.

[Learn More](#)

Parameter	Value
Product Name	VSIG1 Rabbit pAb
Remarks/Alias	GPA34; dJ889N15.1; 1700062D20Rik; VSIG1
Species	Human
GeneID	340547
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-230 of human VSIG1 (NP_872413.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86XK7
Protein Weight	42kDa
Shipping	Ice bag