

Anti-Vsnl1/Vilip-1 Rabbit Monoclonal Antibody For Wb, Elisa - P62760

Item Number: CM0006468

Introduction

Recombinant rabbit monoclonal antibody against human VSNL1/VILIP-1 (Visinin-like protein 1), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying calcium-dependent signaling, neurodegeneration, and Alzheimer's disease biomarkers.

[Learn More](#)

Product Name	VSNL1/VILIP-1 Rabbit mAb
Remarks/Alias	HLP3, VILIP, HPCAL3, HUVISL1, VILIP-1, VSNL1/VILIP-1
Species	Human
GeneID	7447
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2-191 of human VSNL1/VILIP-1 (NP_003376.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62760
Protein Weight	22kDa
Shipping	Ice bag