

Anti-Vsnl1/Vilip-1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P62760

Item Number: CM0007595

Introduction

Rabbit monoclonal antibody against human VSNL1/VILIP-1, validated for WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Detects a band at 22 kDa. Useful in neuroscience and Alzheimer's disease research.

[Learn More](#)

Parameter	Value
Product Name	VSNL1/VILIP-1 Rabbit mAb
Remarks/Alias	HLP3; VILIP; HPCAL3; HUVISL1; VILIP-1; VSNL1/VILIP-1
Species	Human
GeneID	7447
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2-191 of human VSNL1/VILIP-1 (NP_003376.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62760
Protein Weight	22kDa
Shipping	Ice bag