


Anti-Vti1B Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Ueu0

Item Number: CM0027521



Introduction

Rabbit polyclonal antibody against human VTI1B (VTI1L), a v-SNARE involved in endosomal trafficking to lysosomes. Suitable for WB, IHC-P, ELISA; cross-reacts with mouse. Ideal for membrane fusion and lysosomal delivery studies.

[Learn More](#)

Parameter	Value
Product Name	VTI1B Rabbit pAb
Remarks/Alias	VTI1; VTI2; VTI1L; v-SNARE; vti1-rp1; VTI1-LIKE; VTI1B
Species	Human
GeneID	10490
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human VTI1B (NP_006361.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UEU0
Protein Weight	27kDa
Shipping	Ice bag