

Anti-Vwc2L Polyclonal Antibody For Wb And Elisa - Q505H4

Item Number: CM0023001

Introduction

Rabbit polyclonal antibody targeting mouse VWC2L (Brorin-like), validated for Western blot (WB) and ELISA. Cross-reacts with human and mouse. Useful for studies on neurogenesis and bone mineralization. Predicted molecular weight 25 kDa.

[Learn More](#)

Parameter	Value
Product Name	VWC2L Rabbit pAb
Remarks/Alias	Brl; A830006F12Rik; VWC2L
Species	Mouse
Gene ID (Human)	402117
Gene ID	320460
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-50 of mouse VWC2L (NP_796138.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q505H4
Protein Weight	25kDa
Shipping	Ice bag