

Anti-Wapl Rabbit Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - Q7Z5K2

Item Number: CM0009092



Introduction

Rabbit monoclonal antibody against human WAPL, validated for WB, IF/ICC, IP, and ELISA. Reacts with human and mouse. Ideal for studying sister chromatid cohesion and mitosis.

[Learn More](#)

Parameter	Value
Product Name	WAPL Rabbit mAb
Remarks/Alias	FOE, WAPAL, KIAA0261, WAPL
Species	Human
GeneID (Human)	23063
GeneID	23063
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human WAPL (Q7Z5K2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q7Z5K2
Protein Weight	133kDa
Shipping	Ice bag