


Anti-Wasp Rabbit Polyclonal Antibody For Wb, Elisa - P42768

Item Number: CM0016949



Introduction

High-quality rabbit anti-WASP polyclonal antibody (P42768) for Western blot and ELISA applications. Validated in human samples; targets actin nucleation-promoting factor WAS for cytoskeleton and immune research.

[Learn More](#)

Parameter	Value
Product Name	WASP Rabbit pAb
Remarks/Alias	THC, IMD2, SCNX, THC1, WASP, WASPA
Species	Human
GeneID (Human)	7454
GeneID	7454
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-250 of human WASP (NP_000368.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P42768
Protein Weight	53kDa
Shipping	Ice bag