

## Anti-Wave2/Wasf2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Y6W5

Item Number: CM0003195



### Introduction

High-quality rabbit monoclonal antibody specific to human WAVE2/WASF2 (Wiskott-Aldrich syndrome protein family member 2), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for actin dynamics and lamellipodia studies.

[Learn More](#)

| Parameter        | Value                                                                                                                        |
|------------------|------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | WAVE2/WASF2 Rabbit mAb                                                                                                       |
| Remarks/Alias    | IMD2; SCAR2; WASF4; WAVE2; dj393P12.2; WAVE2/WASF2                                                                           |
| Species          | Human                                                                                                                        |
| Gene ID (Human)  | 10163                                                                                                                        |
| Gene ID          | 10163                                                                                                                        |
| Immunogen        | Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 399-498 of human WAVE2/WASF2 (Q9Y6W5). |
| Source           | Rabbit                                                                                                                       |
| Category         | Monoclonal Antibodies                                                                                                        |
| Application      | WB, IHC-P, ELISA                                                                                                             |
| Cross Reactivity | Human, Mouse, Rat                                                                                                            |
| SWISS            | Q9Y6W5                                                                                                                       |
| Protein Weight   | 54kDa                                                                                                                        |
| Shipping         | Ice bag                                                                                                                      |