


Anti-Wbp11 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y2W2

Item Number: CM0018754



Introduction

Rabbit polyclonal antibody targeting human WBP11 (WW domain-binding protein 11) for Western blot, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for splicing mechanism studies.

[Learn More](#)

Parameter	Value
Product Name	WBP11 Rabbit pAb
Remarks/Alias	BUG13; FAP90; NPWBP; SIPP1; VCTRL; CFAP90; VCTERL; WBP-11; PPP1R165; WBP11
Species	Human
Gene ID (Human)	51729
Gene ID	51729
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2W2
Protein Weight	70kDa
Shipping	Ice bag