

Anti-Wbp2NI Rabbit Polyclonal Antibody For Wb - Q6lcg8

Item Number: CM0022564

Introduction

Anti-WBP2NL rabbit polyclonal antibody for Western blot, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Detects 32 kDa WBP2NL (PAWP), involved in fertilization and meiosis. Q6ICG8.

[Learn More](#)

Parameter	Value
Product Name	WBP2NL Rabbit pAb
Remarks/Alias	PAWP; GRAMD7; WBP2NL
Species	Human
GeneID (Human)	164684
GeneID	164684
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human WBP2NL (NP_689826.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q6ICG8
Protein Weight	32kDa
Shipping	Ice bag