


Anti-Wbp4 Rabbit Polyclonal Antibody For Wb, Elisa - 075554

Item Number: CM0021666



Introduction

Rabbit polyclonal antibody against human WBP4 (FBP21), validated for Western blot and ELISA. Cross-reacts with rat. Ideal for studying pre-mRNA splicing and spliceosome function.

[Learn More](#)

Parameter	Value
Product Name	WBP4 Rabbit pAb
Remarks/Alias	FBP21; WBP4
Species	Human
Gene ID (Human)	11193
Gene ID	11193
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-376 of human WBP4 (NP_009118.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	O75554
Protein Weight	43kDa
Shipping	Ice bag