


Anti-Wbscr22 Polyclonal Antibody For Wb, If/Icc, Elisa - O43709

Item Number: CM0030505



Introduction

Rabbit polyclonal antibody targeting human WBSR22 (rRNA methyltransferase), validated for WB, IF/ICC, ELISA. Detects human and rat WBSR22 (~32kDa). Suitable for ribosome biogenesis and chromatin regulation studies.

[Learn More](#)

Parameter	Value
Product Name	WBSR22 Rabbit pAb
Remarks/Alias	WBMT; MERM1; PP3381; HUSSY-3; WBSR22; HASJ4442
Species	Human
GeneID (Human)	114049
GeneID	114049
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1-281 of human WBSR22 (NP_059998.2))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	O43709
Protein Weight	32kDa
Shipping	Ice bag