

Anti-Wbscr27 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8N6F8

Item Number: CM0023976

Introduction

Rabbit polyclonal antibody targeting WBSCR27 (METTL27). Suitable for WB, IHC-P, ELISA with reactivity across human, mouse, and rat. Investigate roles in Williams-Beuren syndrome and methylation pathways.

[Learn More](#)

Parameter	Value
Product Name	WBSCR27 Rabbit pAb
Remarks/Alias	WBSCR27
Species	Human
GeneID (Human)	155368
GeneID	155368
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-245 of human WBSCR27 (NP_689772.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N6F8
Protein Weight	27kDa
Shipping	Ice bag