

Anti-Wdfy4 Polyclonal Antibody For Wb, Elisa - Q6Zs81

Item Number: CM0028235

Introduction

Rabbit polyclonal antibody against human WDFY4, suitable for WB and ELISA. Targets WD repeat- and FYVE domain-containing protein 4, critical for dendritic cell cross-presentation and B-cell autophagy. Validated for human samples.

[Learn More](#)

Parameter	Value
Product Name	WDFY4 Rabbit pAb
Remarks/Alias	C10orf64; WDFY4
Species	Human
GeneID (Human)	57705
GeneID	57705
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human WDFY4 (NP_065996.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q6ZS81
Protein Weight	354kDa
Shipping	Ice bag