


Anti-Wdhd1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - 075717

Item Number: CM0013472



Introduction

Rabbit polyclonal anti-WDHD1 antibody validated for WB, IHC-P, IF/ICC, and ELISA. Recognizes human and rat WDHD1 (126 kDa), a core replisome factor ideal for DNA replication research.

[Learn More](#)

Parameter	Value
Product Name	WDHD1 Rabbit pAb
Remarks/Alias	AND1; CTF4; AND-1; CHTF4; WDHD1
Species	Human
GeneID (Human)	11169
GeneID	11169
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	075717
Protein Weight	126kDa
Shipping	Ice bag