

Anti-Wdr1 Polyclonal Antibody For Wb And Elisa - 075083

Item Number: CM0015758



Introduction

Rabbit polyclonal antibody against human WDR1, validated for WB and ELISA applications. Detects endogenous WDR1 (66 kDa) in human, mouse, and rat samples. Ideal for studying actin cytoskeleton dynamics and cell migration.

[Learn More](#)

Parameter	Value
Product Name	WDR1 Rabbit pAb
Remarks/Alias	AIP1; PFITS; NORI-1; HEL-S-52; WDR1
Species	Human
Gene ID (Human)	9948
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 300-500 of human WDR1 (NP_059830.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075083
Protein Weight	66kDa
Shipping	Ice bag