


Anti-Wdr11 Polyclonal Antibody For Wb, Elisa - Q9Bzh6

Item Number: CM0018389



Introduction

Rabbit anti-WDR11 polyclonal antibody validated for Western blot and ELISA. Targets human WD repeat-containing protein 11, involved in Hedgehog signaling and ciliogenesis. Suitable for research in developmental biology and disease mechanisms.

[Learn More](#)

Parameter	Value
Product Name	WDR11 Rabbit pAb
Remarks/Alias	DR11, HH14, SRI1, BRWD2, WDR15, WDR11
Species	Human
Gene ID (Human)	55717
Gene ID	55717
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 975-1224 of human WDR11 (NP_060587.8)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9BZH6
Protein Weight	137kDa
Shipping	Ice bag