


Anti-Wdr12 Rabbit Polyclonal Antibody For Wb And Elisa - Q9GZL7

Item Number: CM0017996



Introduction

High-quality rabbit polyclonal antibody against human WDR12 (YTM1), validated for Western blot and ELISA. Ideal for ribosome biogenesis research; detects a ~48 kDa band.

[Learn More](#)

Parameter	Value
Product Name	WDR12 Rabbit pAb
Remarks/Alias	YTM1; WDR12
Species	Human
Gene ID (Human)	55759
Gene ID	55759
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human WDR12 (NP_060726.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9GZL7
Protein Weight	48kDa
Shipping	Ice bag