

## Anti-Wdr46 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O15213

Item Number: CM0021692



### Introduction

WDR46 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with Human, Mouse, Rat. Targets nucleolar scaffold protein involved in ribosomal biogenesis. UniProt: O15213.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | WDR46 Rabbit pAb                                                                                                   |
| Remarks/Alias    | UTP7; BING4; FP221; C6orf11; WDR46                                                                                 |
| Species          | Human                                                                                                              |
| Gene ID          | 9277                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 421-610 of human WDR46 (NP_005443.3) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                           |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | O15213                                                                                                             |
| Protein Weight   | 68kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |