

## Anti-Wdr5 Rabbit Polyclonal Antibody For Wb And Elisa - P61964

Item Number: CM0016985



### Introduction

High-quality rabbit polyclonal antibody targeting human WDR5 (WD repeat-containing protein 5), suitable for Western blotting and ELISA. Ideal for chromatin regulation and epigenetic research.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | WDR5 Rabbit pAb                                                                                                  |
| Remarks/Alias    | BIG3; SWD3; BIG-3; CFAP89; WDR5                                                                                  |
| Species          | Human                                                                                                            |
| GeneID (Human)   | 11091                                                                                                            |
| GeneID           | 11091                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-334 of human WDR5 (NP_438172.1). |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Human                                                                                                            |
| SWISS            | P61964                                                                                                           |
| Protein Weight   | 37kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |