

Anti-Wdr6 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nnw5

Item Number: CM0025461

Introduction

Rabbit polyclonal antibody against human WDR6, validated for WB, IF/ICC, and ELISA. Detects endogenous WDR6 in human and mouse samples; involved in tRNA methylation and autophagy regulation.

[Learn More](#)

Parameter	Value
Product Name	WDR6 Rabbit pAb
Remarks/Alias	Trm734; WDR6
Species	Human
Gene ID (Human)	11180
Gene ID	11180
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 800-1000 of human WDR6 (NP_060501.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NNW5
Protein Weight	122kDa
Shipping	Ice bag