


Anti-Wdr61 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Gzs3

Item Number: CM0018129



Introduction

Rabbit polyclonal antibody against human WDR61 (SKIC8/REC14) for WB, IHC-P, ELISA. Cross-reacts with mouse, rat. Targets the 34 kDa PAF1/SKI complex subunit involved in transcriptional regulation and RNA decay. Suitable for epigenetics and cancer research.

[Learn More](#)

Parameter	Value
Product Name	WDR61 Rabbit pAb
Remarks/Alias	SKI8; REC14; WDR61
Species	Human
Gene ID	80349
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-305 of human WDR61 (NP_079510.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZS3
Protein Weight	34kDa
Shipping	Ice bag