


Anti-Wdr77 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Bqa1

Item Number: CM0025958



Introduction

Rabbit polyclonal antibody against human WDR77, validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for methylosome complex and arginine methylation research.

[Learn More](#)

Parameter	Value
Product Name	WDR77 Rabbit pAb
Remarks/Alias	p44; MEP50; MEP-50; HKMT1069; Nbla10071; p44/Mep50; WDR77
Species	Human
GeneID (human)	79084
GeneID	79084
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 153-243 of human WDR77 (NP_077007.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BQA1
Protein Weight	37kDa
Shipping	Ice bag