

Anti-Wdr81 Polyclonal Antibody For Wb, Elisa - Q562E7

Item Number: CM0013634

Introduction

Rabbit polyclonal antibody targeting WDR81, a negative regulator of PI3K activity on endosomal membranes. Validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Supports research into endosomal trafficking and autophagy.

[Learn More](#)

Parameter	Value
Product Name	WDR81 Rabbit pAb
Remarks/Alias	HYC3; CHMRQ; CAMRQ2; SORF-2; PPP1R166; WDR81
Species	Human
GeneID (Human)	124997
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 340-620 of human WDR81 (NP_001157281.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q562E7
Protein Weight	212kDa
Shipping	Ice bag