

Anti-Wee1 Monoclonal Antibody For Ihc-P, Elisa - P30291

Item Number: CM0007328



Introduction

Rabbit monoclonal antibody against human WEE1, suitable for IHC-P and ELISA with mouse cross-reactivity. Ideal for studying G2/M checkpoint regulation and cell cycle control.

[Learn More](#)

Parameter	Value
Product Name	WEE1 Rabbit mAb
Remarks/Alias	WEE1A, WEE1hu, WEE1
Species	Human
Gene ID (Human)	7465
Gene ID	7465
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human WEE1 (NP_003381.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	P30291
Protein Weight	72kDa
Shipping	Ice bag