

Anti-Wfdc1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Hc57

Item Number: CM0021306

Introduction

Rabbit polyclonal antibody targeting WFDC1/PS20, a secreted growth inhibitor. Suitable for WB, IF/ICC, ELISA in human and mouse samples. Ideal for cancer research and protease inhibitor studies.

[Learn More](#)

Parameter	Value
Product Name	WFDC1 Rabbit pAb
Remarks/Alias	PS20; WFDC1
Species	Human
GenelD (Human)	58189
GenelD	58189
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 111-220 of human WFDC1 (NP_067020.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9HC57
Protein Weight	24kDa
Shipping	Ice bag