

Anti-Wfs1 Monoclonal Antibody For Wb, If-P, Elisa - O76024

Item Number: CM0006356

Introduction

This rabbit monoclonal antibody detects human WFS1 (Wolframin) and cross-reacts with rat. It is suitable for WB, IF-P, and ELISA. Useful for studies on ER stress, calcium signaling, and Wolfram syndrome.

[Learn More](#)

Parameter	Value
Product Name	WFS1 Rabbit mAb
Remarks/Alias	WFS; WFRS; WFSL; CTRCT41; WFS1
Species	Human
Gene ID (Human)	7466
Gene ID	7466
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 70-170 of human WFS1 (NP_001139325.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Rat
SWISS	O76024
Protein Weight	100 kDa
Shipping	Ice bag