

Anti-Wfs1 Polyclonal Antibody For Wb, If/Icc, Elisa - 076024

Item Number: CM0023801

Introduction

Rabbit polyclonal antibody against human WFS1 (Wolframin) validated for WB, IF/ICC, and ELISA. Detects human, mouse, and rat WFS1. Ideal for ER stress, calcium homeostasis, and Wolfram syndrome research.

[Learn More](#)

Parameter	Value
Product Name	WFS1 Rabbit pAb
Remarks/Alias	WFS; WFRS; WFSL; CTRCT41; WFS1
Species	Human
GeneID (Human)	7466
GeneID	7466
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-285 of human WFS1 (NP_005996.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human,Mouse,Rat
SWISS	O76024
Protein Weight	100kDa
Shipping	Ice bag