

Anti-Whrn Rabbit Polyclonal Antibody For Wb, Elisa - Q9P202

Item Number: CM0015537



Introduction

Rabbit polyclonal antibody targeting human WHRN (whirlin), applicable for WB and ELISA. Recognizes human and mouse whirlin protein. Suitable for studying hearing, vision, and Usher syndrome.

[Learn More](#)

Parameter	Value
Product Name	WHRN Rabbit pAb
Remarks/Alias	WI; CIP98; USH2D; DFNB31; PDZD7B; WHRN
Species	Human
Gene ID (Human)	25861
Gene ID	25861
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 708-907 of human WHRN (NP_056219.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9P202
Protein Weight	97kDa
Shipping	Ice bag