


Anti-Wif1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y5W5

Item Number: CM0013468



Introduction

WIF1 Rabbit pAb (Wnt inhibitory factor 1) is a rabbit polyclonal antibody validated for Western blot, IHC-P, and ELISA. It reacts with human, mouse, and rat WIF1, a secreted inhibitor of Wnt signaling involved in mesoderm segmentation.

[Learn More](#)

Parameter	Value
Product Name	WIF1 Rabbit pAb
Remarks / Alias	WIF-1; WIF1
Species	Human
Gene ID (Human)	11197
Gene ID	11197
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-220 of human WIF1 (NP_009122.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5W5
Protein Weight	42kDa
Shipping	Ice bag