


Anti-Wif1 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Y5W5

Item Number: CM0010333



Introduction

Rabbit monoclonal antibody targeting human WIF1, validated for WB and ELISA. Recognizes human, mouse, rat WIF1. Ideal for Wnt signaling inhibition studies. Immunogen: synthetic peptide within amino acids 280-379 of human WIF1 (UniProt Q9Y5W5).

[Learn More](#)

Parameter	Value
Product Name	WIF1 Rabbit mAb
Remarks/Alias	WIF-1; WIF1
Species	Human
GeneID (Human)	11197
GeneID	11197
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 280-379 of human WIF1 (Q9Y5W5).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5W5
Protein Weight	42kDa
Shipping	Ice bag