


Anti-Wif1 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - Q9Y5W5

Item Number: CM0007194



Introduction

Rabbit monoclonal antibody against human WIF1, validated for WB, IF-P, and ELISA. Cross-reacts with mouse. Targets Wnt inhibitory factor 1, a secreted WNT signaling antagonist. Ideal for developmental and cancer research.

[Learn More](#)

Parameter	Value
Product Name	WIF1 Rabbit mAb
Remarks/Alias	WIF-1; WIF1
Species	Human
Gene ID	11197
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-379 of human WIF1 (NP_009122.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y5W5
Protein Weight	42kDa
Shipping	Ice bag